



Never adjust inside of the crown (use fit-checker “black”)

### Cementation Science

	Zirconium Oxide	Lithium disilicate glass-ceramic (Emax, Celtra and etc.)
<b>Etching/Cleaning</b> (inside crown)	Ivoclean *** (rinse out after 20 sec.and dry)	5% HF Acid or Ivoclean *** (rinse out after 20 sec. and dry)
<b>Conditioning/ Silanization</b> (inside crown)	Monobond® Plus (blow thin after 60 sec.)	Monobond® Plus (blow thin after 60 sec.)
<b>Cementation System</b> (on preparation)	Multilink® Automix Unicem	Multilink® Automix SpeedCEM Viva Glass® CEM

\*\*\*See Ivoclean on Youtube

**Modern Dental Technology Laboratory** • 1295 Leila Avenue, Winnipeg, MB R2P 2Z3, Canada Phone  
number: +1 (204) 632-0407 • cell: +1 (204) 805-2589  
email address: [Moderdentallab@hotmail.com](mailto:Moderdentallab@hotmail.com)